

RE-C 6 0 0 Controller Specifications for Ring Pump series RP-Q-B / RP-Q-C / RP-TX / RP-HX Bi-polar Stepper Motors


Product Name	Low voltage bi-polar stepper motor controller
Model Number	RE-C600 : controller unit
Applicable Pumps	Ring Pump RP-Q-B series, RP-Q-C series Ring Pump RP-TX series Ring Pump RP-HX series
Drive Voltage	RP-Q-B / RP-Q-C Series: VM=3.0 V RP-TX / RP-HX Series : VM=2.6 V
Input Voltage	RE-C600 : DC5V, 1.0 A (AC-DC adaptor of AC100-240V 50/60Hz is included as accessory)
Control Method	Voltage drive method for bi-polar coil
Excitation Method	1/4 microstep drive method
Motor Speed	1 - 2,000 pps (4 - 8,000 Hz) Digital frequency division
Functions	Various controls using PC software SINGL-MODE Program-Mode

Connection Terminals: RE-C600 (Size: 100 x 100 x 34 mm)



- USB Connection Terminal (USB A-microB: Included as accessory)
- Pump Connection Terminals-1
- Pump Connection Terminals-2
- DC-IN Connection Terminal (AC Adaptor: Included as accessory)
- Start Switch
- Start-LED (BLUE)
- Pause-Switch
- Pause-LED (GREEN)
- Power-LED (RED)

Pump wiring to motor connection terminal



- Accessories Included: RE-C600 :**
- Type A-micro B USB Cable
 - USB Memory Stick (PC Software)
 - AC-DC Adapter

Motor Speed (PPS) and Motor Rotation for Each Series

RP-Q-B Series (1/150.95 gear ratio motor specification)

Motor Speed PPS	1	5	10	50	100	200	300	400	500	800	1000	1200	1250
Motor Rotation rpm	0.02	0.10	0.20	0.99	1.99	3.98	5.96	7.96	9.94	15.92	19.87	23.84	24.83

7.96

RP-Q-C Series (1/51.45 gear ratio motor specification)

Motor Speed PPS	1	5	10	50	100	200	300	400	500	800	1000	1100	
Motor Rotation rpm	0.06	0.29	0.58	2.92	5.83	11.66	17.49	23.32	29.15	46.64	58.31	64.14	

RP-TX Series (1/135.8 gear ratio motor specification)

Motor Speed PPS	1	5	10	50	100	200	300	400	500	800	1000	1200	1500
Motor Rotation rpm	0.02	0.11	0.22	1.10	2.21	4.42	6.63	8.84	11.05	17.67	22.09	26.51	33.14


RP-HX Series (1/50 gear ratio motor specification)

Motor Speed PPS	1	5	10	50	100	200	300	400	500	800	1000		
Motor Rotation rpm	0.06	0.30	0.60	3.00	6.00	12.00	18.00	24.00	30.00	48.00	60.00		

Special Notes:

- a) Discharge Rate: Please input "Discharge Rate per Rotation ($\mu\text{L}/\text{rev}$)".
(Refer to data included with the pumps)
- b) Control the motors within the range to prevent motor stepping-out (or motor stopping or vibrating).

SINGLE



<Single Settings>

- ① Pump Selection
 - Click the "v" (drop down) and select the pump.
 - RP-Q-B RP-Q-C
 - RP-TX RP-HX
 - Operating Voltage: VM value is selected automatically.
 - RP-Q-B/RP-Q-C VM = 3.0V (fixed)
 - RP-TX/RP-HX VM = 2.6V (fixed)
- ② DIRECTION
 - rotating direction. (Clockwise CW or Counter-Clockwise CCW)
- Pump Speed
 - 1 ~ 2000 PPS (Limit: 1500)
- FREE
 - Click FREE to supply media (solution) for pump
 - Click STOP.

<PROGRAM Settings>

- Input Resolution "μL/rev"
 - Input the value from "Discharge Volume per Rotation (μL/rev)" that comes attached with the pump.
- ENTER
 - Click ENTER to fix the "Resolution" settings.
- ROTATION
 - 1 ~ 30
- Flow Rate (uL)
 - $\text{Flow Rate (uL)} = \text{Resolution} \times \text{Rotation}$
- START
 - START button while connected to the PC (computer) with USB.
 - * Auto Stop
- EXIT
 - Disconnect from PC (computer).

PROGRAM

<Initial Settings>

- Pump Selection Click the "v" (drop down) and select the pump.
RP-Q-B RP-Q-C
RP-TX RP-HX
- VOLTAGE Operating Voltage: VM value is selected automatically.
RP-Q-B/RP-Q-C VM = 3.0V (fixed)
RP-TX/RP-HX VM = 2.6V (fixed)
- Fast Forward Operation - 1 Click the "v" (drop down) and rotating direction. (Clockwise CW or Counter-Clockwise CCW)
- Fast Forward Operation - 2 START/STOP the optimal high speed pump operation to supply the media (solution).
* Click START to supply media (solution) for each pump on each channel. Once the media (solution) flows to the tip of the discharge tube, click STOP.

<PROGRAM Settings>

- Input Resolution "μL/rev" Input the value from "Discharge Volume per Rotation (μL/rev)" that comes attached with the pump.
- Specify Flow Rate "μL/min" Input desired flow rate (μL/min).
- ENTER Click ENTER to fix the "Resolution" and "Flow Rate" settings. The optimal "PPS" and "RPM" values will be automatically displayed.
- Select DIRECTION Select between Clockwise (CW) or Counter-Clockwise (CCW)
- Operation Time Set operation (ON) time
* Setting range between 00h 00m 01s to 99h 59m 59s
- Stop Time Set stop (OFF) time
* Setting range between 00h 00m 01s to 99h 59m 59s
(3)
- POS-CLEAR Select between POS 1 to POS 5 () and click POS-CLEAR to clear setting for that POS. (To reset settings)
- RESET RESET all settings in POS 1 to POS 5.
- Specify POS Click the "v" (drop down) and select POS.
Allows up to 5 different program condition settings in POS 1 to POS 5. Program is reflected on the POS by repeating to .
(Flow Rate, ENTER, DIRECTION, Operation Time, Stop Time)
- REPEAT Input number of repeat operation for specified program table (POS 1 to POS 5).
Initial setting is "1". Can input value between 1 to 99.
- WAIT TIME (h:m) If needed, specify the wait time before the START time for each channel. (Input h:m)

<WRITE Settings>

- WRITE Click WRITE to input the PROGRAM settings to controller's ROM (memory) (During WRITE, it will be in WAIT mode)

<Operation Settings>

"Controller Switch - START"

Use the controller's switch to START / STOP / PAUSE without the PC.

"PC-START"

START

START / STOP button while connected to the PC (computer) with USB.

* Each click will change between START STOP.

PAUSE

PAUSE button while connected to the PC (computer) with USB.


* Each click will change between PAUSE RESTART PAUSE.

EXIT

Disconnect from PC (computer).

NB. When controlling from PC side, "START" and "Pause" of main unit SW must be set to <OFF> position to prevent possible unstable operation.

SAVE-FILE / OPEN FILE



SAVE-FILE

Click on SAVE-FILE to save the PROGRAM setting (Type in file name)

OPEN-FILE

Click on OPEN-FILE to open saved PROGRAM setting.

* Once the PROGRAM setting is opened, click WRITE to save it in the computer ROM.

* Click all the CH that needs to operate (letters turn from grey to black), and click CH-SET to save it in the computer ROM (memory).